

Oxford Nanopore sequencing of pharmacogenomic (PGx) loci delivers unambiguous variant calling in a single assay

Featuring highly complex structural variants (SVs), pseudogene homology, and copy number variants (CNVs), PGx loci are incredibly challenging to resolve using conventional microarrays or short-read sequencing. Oxford Nanopore sequencing provides the solution to resolve these complexities, by providing any-length reads for calling all variants in a single assay, via an accessible, scalable platform.

The following highlighted publications showcase how researchers are taking advantage of this technology, focusing on targeted analysis through either sample prep-based enrichment, or adaptive sampling — a fast and flexible method unique to Oxford Nanopore devices that involves software-controlled enrichment.

Sample prep-based targeted enrichment

Featured publication	Workflow	Research summary
<p>Characterisation of complex structural variation in the <i>CYP2D6-CYP2D7-CYP2D8</i> gene loci using single-molecule long-read sequencing</p> <p>Turner, A.J. et al. <i>Front. Pharmacol.</i> 14:1195778 (2023).</p>	<p>Enrichment method: CRISPR-Cas9 (PCR-free)</p> <p>Device: MinION™</p>	<p>Aim: develop an improved method to characterise complex SVs and CNVs in the <i>CYP2D6-CYP2D7-CYP2D8</i> loci, in a single assay.</p> <p>Results: the researchers achieved full haplotype phasing and precise characterisation of <i>CYP2D</i> SVs. Novel <i>CYP2D6</i> suballeles were identified, and complex diplotypes were resolved in multiple samples. The team uncovered 17 <i>CYP2D7</i> and 18 <i>CYP2D8</i> unique haplotypes. They also determined the haplotype structure of the entire region for all samples.</p>
<p>Approach for phased sequence-based genotyping of the critical pharmacogene dihydropyrimidine dehydrogenase (<i>DPYD</i>)</p> <p>Ambrodjii, A. et al. <i>Int. J. Mol. Sci.</i> 25(14):7599 (2024).</p>	<p>Enrichment method: PCR</p> <p>Device: GridION™</p>	<p>Aim: develop a cost-effective method to genotype the complex <i>DPYD</i> gene.</p> <p>Results: all <i>DPYD</i> variants detected <i>a priori</i> were identified. Fully phased genotypes were obtained for all samples, overcoming the limitations of conventional genotyping methods. An additional 11 heterozygous genotypes — previously missed by standard methods — were detected, then validated by orthogonal methods. As few as 619 high-quality, full-length reads were needed to reliably call and phase variants.</p>

Adaptive sampling-based targeted enrichment

Featured publication	Workflow	Research summary
<p>Targeted haplotyping in pharmacogenomics using Oxford Nanopore Technologies' adaptive sampling</p> <p>Deserranno, K. et al. <i>Front. Pharmacol.</i> 14:1286764 (2023).</p>	<p>Enrichment method: Adaptive sampling</p> <p>Device: PromethION™ 24</p>	<p>Aim: assess the potential of adaptive sampling on an extended panel of 1,036 genes associated with drug response, to determine its feasibility for cost-effective PGx gene analysis.</p> <p>Results: calling and phasing of single nucleotide variants (SNVs), indels, and SVs was performed successfully. Even with three samples multiplexed on a single flow cell, high recall (99.35%) and precision (99.84%) were obtained for calling of targeted PGx variants.</p>
<p>Cas9 targeted nanopore sequencing with enhanced variant calling improves <i>CYP2D6-CYP2D7</i> hybrid allele genotyping</p> <p>Rubben, K. et al. <i>PLoS Genet.</i> 18(9):e1010176 (2022).</p>	<p>Enrichment method: Adaptive sampling + CRISPR-Cas9</p> <p>Device: GridION</p>	<p>Aim: develop a more accurate, complete method, versus current approaches, for genotyping complex PGx regions such as the highly polymorphic <i>CYP2D6</i> gene.</p> <p>Results: <i>CYP2D6-CYP2D7</i> hybrid alleles were accurately identified via an optimised Oxford Nanopore sequencing-based assay and tailored analysis pipeline for simultaneous SNV and SV detection. Novel SVs and SNVs were also detected that were missed by other genotyping methods.</p>

'We conclude that PGx based on targeted [long-read sequencing] is a valuable tool to advance the implementation of personalized medicine.' Deserrano *et al.* 2023






Find out more about pharmacogenomic sequencing with Oxford Nanopore: nanoporetech.com/pgx

Additional references

1. Anukul, N. et al. Ultrarapid and high-resolution HLA class I typing using transposase-based nanopore sequencing applied in pharmacogenetic testing. *Front. Genet.* 14:1213457 (2023). DOI: <https://doi.org/10.3389/fgene.2023.1213457>
2. Hitchman, L.M. et al. Allelic diversity of the pharmacogene *CYP2D6* in New Zealand Māori and Pacific peoples. *Front. Genet.* 13:1016416 (2022). DOI: <https://doi.org/10.3389/fgene.2022.1016416>
3. Payen, L. Presentation. Available at: <https://nanoporetech.com/resource-centre/importance-of-adaptive-sampling-in-nanopore-long-read-sequencing-in-pharmacogenetics> [Accessed 19 February 2025]
4. Verma, R. et al. A nanopore sequencing-based pharmacogenomic panel to personalise tuberculosis drug dosing. *Am. J. Respir. Crit. Care Med.* 209(12):1486-1496 (2024). DOI: <https://doi.org/10.1164/rccm.202309-1583OC>
5. Zhang, L.-Q. et al. Genetic variants, haplotype determination, and function of novel alleles of *CYP2B6* in a Han Chinese population. *Heliyon.* 10(7):e28952 (2024). DOI: <https://doi.org/10.1016/j.heliyon.2024.e28952>



phone +44 (0)845 034 7900
email support@nanoporetech.com
nanoporetech.com

 [oxford-nanopore-technologies](https://www.linkedin.com/company/oxford-nanopore-technologies)
 [@nanopore](https://twitter.com/nanopore)
 [@nanoporetech.com](https://www.facebook.com/nanoporetech.com)