



Assembling bacterial genomes using long nanopore sequencing reads

In order to understand the true diversity and biology of microorganisms, producing fully annotated, complete genomes is essential. However, 90% of bacterial genomes are predicted to be incomplete1.

Long, PCR-free nanopore sequencing reads enable the assembly of complete, referencequality microbial genome sequences. Nanopore sequencing shows a lack of bias in GC-rich regions, in contrast to other sequencing platforms², and can span repeat-rich sequences and structural variants that are inaccessible to traditional sequencing technologies.

Here we present a simple workflow for bacterial genome assembly from a single-organism culture, using MinION™ Flow Cells on MinION or GridION™ sequencing devices.




EXTRACTION:

obtaining high molecular-weight DNA


To extract DNA from a bacterial isolate, we recommend the QIAGEN Genomic-tip 500/G. which we have found to produce the longest read lengths and highest yield. Performing size selection on the extracted DNA further increases read length N50; size selection options include Agencourt AMPure XP beads and the Oxford Nanopore Short Fragment Eliminator Expansion.

Find more extraction protocols, such as those for stool, soil, and culture, as well as alternative size selection methods: community.nanoporetech.com/docs/prepare



LIBRARY PREPARATION:

selecting a kit



To prepare gDNA for sequencing, we recommend the Ligation Sequencing Kit, providing the greatest throughput and control over read lengths.

Find out more about library preparation kits, including rapid,

store.nanoporetech.com/sample-prep.html

Multiplexing options are available to increase the cost efficiency of your sequencing. We recommend the Native Barcoding Kits, which are PCR-free, and enable up to 96 samples to be sequenced in multiplex on one flow cell. Alternatively, PCR-based options are also available for the multiplexed sequencing of up to 96 samples.




SEQUENCING: achieving ultimate flexibility with MinION Flow Cells



We recommend sequencing bacterial genomes on MinION Flow Cells. These flow cells can be used individually on the portable MinION and MinION Mk1C sequencing platforms; alternatively, the benchtop GridION enables on-demand sequencing on up to five MinION Flow Cells. You can therefore adjust flow cell number according to the degree of multiplexing, number of samples, and/or your experimental goals.

Find out more about MinION: nanoporetech.com/products/minion



For assembly, we recommend basecalling in high accuracy mode, and sequencing to a minimum depth of $30 \times 0.00 \times 0.000 \times 0.0000$ kb reads per genome. We suggest multiplexing 12–24 samples per flow cell, scaled up or down according to your samples (e.g. expected genome length; fragment length distribution) and experimental aims. However, if variant calling is also desired, increasing sequencing depth is advisable.

Find out more and compare nanopore sequencing platforms: nanoporetech.com/products/comparison

ANALYSIS:

selecting an assembly tool

To perform bacterial genome assembly, we suggest using the third-party *de novo* assembly tool **Flye**³. This analysis package represents a complete pipeline, taking raw nanopore reads as input, and producing polished contigs as output. We also recommend one round of polishing with **Medaka**⁴. These tools can be found on GitHub.

Bacterial genomes up to ~10 Mb can be assembled with Flye using a standard laptop (~16 GB memory).

This complete analysis pipeline is also available in EPI2ME Labs, which provides best practice workflows and interactive tutorials to support the analysis of your nanopore sequencing data and develop your bioinformatics skills. Find out more at labs.epi2me.io.

Find out more about data analysis solutions: nanoporetech.com/analyse

FASTQ input file Assembly Flye Polishing Medaka (1 round)

FASTA output file

Find out more at: nanoporetech.com/applications/microbiology



Twitter: @nanopore www.nanoporetech.com

References:

- Land, M. et al. Insights from 20 years of bacterial genome sequencing. Funct. Integr. Genomics. 15(2): 141–161 (2015).
- Browne, P. D. et al. GC bias affects genomic and metagenomic reconstructions, underrepresenting GC-poor organisms. GigaScience. 9(2): giaa008 (2020).
- Kolmogorov, M. et al. Assembly of long, error-prone reads using repeat graphs. Nat. Biotech. 37: 540-546 (2019).
- 4. GitHub. Medaka. Available at: https://github.com/nanoporetech/medaka [Accessed: 25 August 2022].

Oxford Nanopore Technologies and the Wheel icon, GridION, and MinION are registered trademarks of Oxford Nanopore Technologies plc in various countries. All other brands and names contained are the property of their respective owners. © 2022 Oxford Nanopore Technologies plc. All rights reserved. Oxford Nanopore Technologies products are not intended for use for health assessment or to diagnose, treat, mitigate, cure, or prevent any disease or condition. WF_1067(EN)_V3_05Sep2022.